



# Amtryke Therapeutic Tricycle Request Form

(Must be filled out completely by adult rider or parent/guardian)

Recipient's Name: \_\_\_\_\_ Age: \_\_\_\_\_ Date of Request: \_\_\_\_\_  
Mailing Address: \_\_\_\_\_ Phone #: \_\_\_\_\_  
City/State/Zip: \_\_\_\_\_ Email: \_\_\_\_\_  
Diagnosis: \_\_\_\_\_

## If Recipient is Under Age 18

Parent/Guardian Name: \_\_\_\_\_  
**If different from above**  
Mailing Address: \_\_\_\_\_ Phone #: \_\_\_\_\_  
City/State/Zip: \_\_\_\_\_ Email: \_\_\_\_\_

Secondary Contact Name: \_\_\_\_\_ Phone #: \_\_\_\_\_

Treating Therapist's Name: \_\_\_\_\_  
Phone #: \_\_\_\_\_ Email: \_\_\_\_\_

How did you hear about the Amtryke Therapeutic Tricycle? (Check all that apply)

Therapist  Website  AMBUCS Member  Other: \_\_\_\_\_

Will you need financial assistance to obtain the tricycle?  Yes  No

If yes, how much can you pay? \_\_\_\_\_

**Note:** Amtryke therapeutic tricycles are distributed based on available funds and need. Individual placements of Amtryke therapeutic tricycles are at the discretion of the local chapter or parent organization.

Tell us about the recipient\*:  
\_\_\_\_\_  
\_\_\_\_\_  
\_\_\_\_\_

\*This information will be made public to help obtain funding. Please don't include information you don't want shared.

Including a photo of the recipient will help us obtain funding to help pay for the Amtryke more quickly. Digital images preferred but we also accept professionally printed glossy photos. No photocopies or folded images.

By including a photo, you are giving consent for AMBUCS to use the image online and in print to help obtain funding.

Signature: \_\_\_\_\_ Date: \_\_\_\_\_

**Individuals will not be considered for placement until all three forms are returned to AMBUCS: this form filled out by the recipient/guardian, as well as the Assessment Form and Tryke Selection Form filled out by the therapist.**

# Amtryke Therapeutic Tricycle Waiver Form

(Must be filled out completely by adult rider or parent/guardian)

AMBUCS members nationwide are dedicated to creating opportunities for mobility and independence for people with disabilities by providing AmTryke therapeutic tricycles, offering educational scholarships to therapy students and performing various forms of community service.

**Purpose:** The AmTryke therapeutic tricycle is designed for people with disabilities. It creates a feeling of freedom, builds self-esteem, strengthens muscles and improves motor coordination and range of motion—all while making exercise fun.

**Steering:** Initially, the rider may have difficulty turning or changing directions. Encourage the rider to go straight ahead, back up and slowly turn around. On many models there are three steering options for the Amtryke. On the front column of the tricycle you will find two holes for the steering pin. The top hole is straight steering, the bottom hole allows a 20 degree turning radius. Leaving the pin out gives the rider free steering.

## Safety Cautions

- Fast speeds and sharp turns can cause the Amtryke therapeutic tricycle to tip or turn over.
- Always wear a helmet when riding an Amtryke. Use of other protective gear is highly recommended.
- Adult supervision required if used by younger or developmentally delayed riders.
- Use caution near vehicles, swimming pools and other bodies of water, hills, alleys and sloped driveways.
- Always wear shoes.
- Never allow more than one rider.
- Use of the steering pin is recommended to prevent over-steering and possible tip-overs.

*The information contained in this document is not intended nor implied by National AMBUCS™, Inc. to be professional medical advice by National AMBUCS, Inc. Always seek the advice of your physician, therapist or other qualified healthcare provider prior to starting any treatment or with any question you may have regarding a medical condition. Nothing contained in this document is intended by National AMBUCS, Inc. to be for medical diagnosis or treatment by National AMBUCS, Inc. or on behalf of National AMBUCS, Inc.*

*In no event shall National AMBUCS, Inc. be liable for any direct, indirect, incidental, consequential, special, exemplary, punitive, or any other monetary or other damages, personal injury or property damages, fees, fines, costs, attorney fees, or liabilities of any kind arising out of or relating in any way to this service or use of the Amtryke® therapeutic tricycle, and/or content or information provided herein.*

I agree that I may be photographed by National AMBUCS. I also agree that my photo and name may be used in promotional efforts for National AMBUCS and Amtryke, LLC. I further grant National AMBUCS the ability to use the photos and name for advertising/publicity purposes without additional compensation, except where prohibited by law.

**By signing below, I acknowledge that I have read and understood this liability waiver.**

Recipient's Name (printed): \_\_\_\_\_

Adult Recipient Signature: \_\_\_\_\_

**If Recipient is Under Age 18**

Legal Guardian Name: \_\_\_\_\_

Legal Guardian Signature: \_\_\_\_\_ Date: \_\_\_\_\_

**Individuals will not be considered for placement until all three forms are returned to AMBUCS: this Form filled out by the recipient/guardian, as well as the Assessment Form and Tryke Selection Form filled out by the therapist.**


# Amtryke Assessment Form

(Must be filled out completely by therapist)

Recipient's Name: \_\_\_\_\_  
 Age: \_\_\_\_\_ Weight (lbs.):\* \_\_\_\_\_ Height (inches):\* \_\_\_\_\_  
 Diagnosis(es):\* \_\_\_\_\_

\*This private information is only used to help appropriately fit the rider.

RIDER'S MEASUREMENTS				
Arm Measurements (inches)				Total Length
Left	A to B:		B to C:	
Right	A to B:		B to C:	
Trunk	A to D:			
Leg Measurements (inches)				Total Length
Left	D to E:		E to F:	
Right	D to E:		E to F:	



A	Center of Shoulder
B	Center of Elbow
C	Center of Digit Crease
D	Center of Hip
E	Center of Knee
F	Bottom of Foot

**Arm Length & Leg Length Measurements are critical to correct Amtryke Selection**

Sizing Chart is available online:  
[www.ambucs.org/riders/wish-list/sizing-chart/](http://www.ambucs.org/riders/wish-list/sizing-chart/)

**Notes on Provided Measurements (if any):**

Helmet Sizing	
Size	Measurement (head circumference)
XXS	18.5" to 19.5"
XS	20.5" to 22"
S/M	22" to 23.6"
L/XL	23.6" to 35.75"

Therapist Name: \_\_\_\_\_ Is this the treating therapist?  Yes  No

Credentials: \_\_\_\_\_

Phone: \_\_\_\_\_ Email: \_\_\_\_\_

Facility Name: \_\_\_\_\_

Street Address: \_\_\_\_\_ City: \_\_\_\_\_ State: \_\_\_\_\_ Zip: \_\_\_\_\_

Is this facility an Amtryke Evaluation Site?  Yes  No  Not Sure

Therapist comments concerning recipient or goals:

This request/assessment is directed to:

- Local AMBUCS Chapter Name: \_\_\_\_\_  
 National Wish List (AMBUCS Resource Center)

By signing below, you are signifying that in your professional opinion this rider would benefit from an Amtryke. You assume no liability.

Therapist Signature: \_\_\_\_\_ Date: \_\_\_\_\_

Ship Amtryke To			
Name/Facility: _____	Phone: _____		
Street Address: _____			
City: _____	State: _____	Zip: _____	

Recipient's Name:

### Thanks for choosing an Amtryke therapeutic tricycle!

In order to accommodate the widest variety of people and disabilities, Amtryke offers many tryke models and each can be customized in a variety of ways. Following the steps below will help you choose to the perfect tryke for your client from what might seem like a dizzying array of options.

Remember you can always refer to our website, [www.amtrykestore.com](http://www.amtrykestore.com), or the Amtryke catalogue for more information and product images.

**Step 1:** Fill out the first page of the Amtryke Assessment Form.

**Step 2:** Choose the way the tryke will be propelled: **Hand & Foot, Foot, or Hand**. Your choice should be based on the rider's ability and therapy goals.

**Hand & Foot** trykes improve coordination, strength and range of motion. Using all four extremities helps with weakness in any area, even general weakness, and can positively affect tone.

**Foot** trykes were developed in response to requests from therapists for a traditional tricycle for riders with special needs.

**Hand** trykes are designed for persons whose lower limbs lack function or those who need special therapy for the upper extremity.

**Step 3:** Take rider's measurements from the front of this Assessment Form and compare them to our Sizing Chart. (This will narrow the choices considerably)

**Step 4:** Choose a drive. (If it doesn't mention a choice, then ignore this step.)

Amtrykes come with two drive possibilities: **fixed** drive or **geared** drive. Tryke models have been carefully designed so the drive matches the therapeutic goals of the equipment; therefore all drives are not available on all trykes.

A **fixed** drive, commonly known as a 'fixie,' works on a mechanical level to help individuals make a full pedal rotation. The foot crank is constantly in motion for full therapeutic effect. Coasting is not possible with a positive drive; when limb motion stops, the bike does as well.

Riders of **geared** drive trykes must be able to make a complete pedal rotation on their own. This tryke is suitable for riders who need help with balance and a stable sitting position and who have the cognitive ability to steer successfully and apply the coaster brake or hand brake. The key feature of geared trykes is the ability to coast.

**Step 5:** Chose any adaptations and/or accessories needed by the rider. *Each tryke model can only be customized in the ways noted in its own model section in the Tryke Selection Forms or with the generic accessories listed below. If a customization option is not listed, it is because of design or other practical constraints in offering it on a particular model.*

### GENERIC ACCESSORIES (not model specific)

<b>Safety Items</b>	<input type="checkbox"/> <b>Safety Flag</b>	<input type="checkbox"/> <b>NEW! Medical Cart</b> (carries O <sub>2</sub> bottle/ventilator/tubing)	<input type="checkbox"/> <b>Water Bottle w Cage</b>	<input type="checkbox"/> <b>License Plate</b>
<b>Leg &amp; Foot Items</b>	<input type="checkbox"/> <b>Foot Cups (pair):</b> <input type="checkbox"/> Small <input type="checkbox"/> Medium	<input type="checkbox"/> <b>Knee Adductor Strap:</b> <input type="checkbox"/> Small <input type="checkbox"/> Large <input type="checkbox"/> Medium	<input type="checkbox"/> <b>Pedal Block (1 = ¾")</b> ____ qty	
<b>Hand Items</b>	<input type="checkbox"/> <b>Variable Range of Motion Kit</b> (only for Hand & Foot Cycles)	<input type="checkbox"/> <b>Wrist Wraps</b> (Includes right & left) <input type="checkbox"/> X-Small <input type="checkbox"/> Large <input type="checkbox"/> Small <input type="checkbox"/> Medium	<input type="checkbox"/> <b>Wrist Brace Mitt:</b> <input type="checkbox"/> <b>Right</b> <input type="checkbox"/> <b>Left</b> <input type="checkbox"/> XX-Small <input type="checkbox"/> Medium <input type="checkbox"/> X-Small <input type="checkbox"/> Small	

Recipient's Name:

## FOOT TRYKES

<input type="checkbox"/> <b>New! 1410 FOOT TRYKE</b> – (with 11" wheels, fixed drive) <b>50-FC-0100</b> Arm Length 7-20", Leg Length 14-21"
<input type="checkbox"/> Option #1 – Standard Seating System: Blue Bucket Seat
<input type="checkbox"/> Option #2 – Upgraded Seating System: Snappy Seat System Snappy Seat Accessories: <input type="checkbox"/> Notched Seat Bottom Cushion* <input type="checkbox"/> Notched Wide Bottom Cushion* <input type="checkbox"/> XL Seat Back <input type="checkbox"/> Laterals <input type="checkbox"/> H-Harness: <input type="checkbox"/> 8" <input type="checkbox"/> 11.5" *Cannot be used with rear steering
1410 Accessories: <input type="checkbox"/> Rear Steering Kit <input type="checkbox"/> Pull/Steering Bar <input type="checkbox"/> Push Bar <input type="checkbox"/> Heavy-Duty Push Bar

<input type="checkbox"/> <b>SNAPPY FOOT TRYKE</b> (with 12" wheels, fixed drive) – <b>50-FC-0150</b> Arm Length 9-13", Leg Length 17-22"
<input type="checkbox"/> Option #1 – Standard Seating System: Snappy Seat System Snappy Seat Accessories: <input type="checkbox"/> Notched Seat Bottom Cushion <input type="checkbox"/> Notched Wide Bottom Cushion <input type="checkbox"/> XL Seat Back <input type="checkbox"/> Laterals <input type="checkbox"/> H-Harness: <input type="checkbox"/> 8" <input type="checkbox"/> 11.5"
<input type="checkbox"/> Option #2 – Alternate Seating System: Blue Bucket Seat
<input type="checkbox"/> Option #3 – Alternate Seating System: Saddle Seat/Gray Plastic Back/Seat Post
Snappy Accessories: <input type="checkbox"/> Pull/Steering Bar <input type="checkbox"/> Push Bar <input type="checkbox"/> Heavy-Duty Push Bar <input type="checkbox"/> Exercise Pedals

<input type="checkbox"/> <b>1412 ProSeries FOOT TRYKE</b> (with 12" wheels, fixed drive) – <b>50-FC-1412</b> Arm Length 14-22", Leg Length 21-26"
<input type="checkbox"/> Option #1 – Standard Seating System: Saddle Seat/1600 Simple Seat Back Seat Bottom Alternates: <input type="checkbox"/> Pommel Saddle Seat <input type="checkbox"/> Medium (11.5x10") <input type="checkbox"/> Skinny Saddle Seat (8.5x6") Trunk Support Upgrades: <input type="checkbox"/> ProSeries Seatback ↴ <input type="checkbox"/> ProSeries Head Rest <input type="checkbox"/> ProSeries Full Padded Back (includes push grip, backpad & 2 laterals) <input type="checkbox"/> ProSeries Lumbar Pad <input type="checkbox"/> ProSeries Back Spacer Kit <input type="checkbox"/> H-Harness: <input type="checkbox"/> 8" <input type="checkbox"/> 11.5" <input type="checkbox"/> AM Full Padded Back
<input type="checkbox"/> Option #2 – Upgraded Seating System: Gray Bucket Seat Bucket Seat Accessories: <input type="checkbox"/> H-Harness: <input type="checkbox"/> 8" <input type="checkbox"/> 11.5"
<input type="checkbox"/> Option #3 – Upgraded Seating System: Snappy Seat System Snappy Seat Accessories: <input type="checkbox"/> Notched Seat Bottom Cushion <input type="checkbox"/> Notched Wide Bottom Cushion <input type="checkbox"/> XL Seat Back <input type="checkbox"/> Laterals <input type="checkbox"/> H-Harness: <input type="checkbox"/> 8" <input type="checkbox"/> 11.5"
1412 Accessories: <input type="checkbox"/> Pull/Steering Bar <input type="checkbox"/> Push Bar <input type="checkbox"/> Heavy-Duty Push Bar <input type="checkbox"/> Knee Separator: <input type="checkbox"/> 2" <input type="checkbox"/> 5" <input type="checkbox"/> Toe Pedal Leveler Pulley <input type="checkbox"/> Hand Brake <input type="checkbox"/> Exercise Pedals <input type="checkbox"/> Hi-Rise Handlebars: <input type="checkbox"/> 5" (BMX) <input type="checkbox"/> 8" <input type="checkbox"/> Rear Steering Kit ( <b>NEW! Takes the place of 1500 Series</b> )





**Recipient's Name:**

**JT-2000 FOOT TRYKE** (Recumbent 21 Speed) – **50-FC-2000** (advanced riders only)

JT-2000 Accessories:       Exercise Pedals       XL Exercise Pedals       Expanding Pedals       Digital Speedometer  
                                  Dual Hand Brake                       Toe Clips                                       Rearview Mirror

**JT-2300-USS FOOT TRYKE** (Recumbent 21 Speed with **Under Seat Steering**) – **50-FC-2300-USS** (advanced riders only)

JT-2300-USS Accessories:       Exercise Pedals       XL Exercise Pedals       Expanding Pedals       Digital Speedometer  
                                  Dual Hand Brake                       Toe Clips                                       Rearview Mirror

Therapist Assembly Notes or Comments (if any):



Recipient's Name:

## HAND-FOOT TRYKES

<input type="checkbox"/> <b>New! AM-10 HAND-FOOT TRYKE</b> (with 4" crank arms) <b>50-HFC-0105</b> Arm length 12-16", leg length 14-21"
<input type="checkbox"/> Option #1 – Standard Seating System: Blue Bucket Seat
<input type="checkbox"/> Option #2 – Upgraded Seating System: Snappy Seat System Snappy Seat Accessories: <input type="checkbox"/> Notched Seat Bottom Cushion* <input type="checkbox"/> Notched Wide Bottom Cushion* <input type="checkbox"/> XL Seat Back <input type="checkbox"/> Laterals <input type="checkbox"/> H-Harness: <input type="checkbox"/> 8" <input type="checkbox"/> 11.5"
AM-10 Accessories: <input type="checkbox"/> Rear Steering Kit <input type="checkbox"/> Push Bar <input type="checkbox"/> Heavy-Duty Push Bar <input type="checkbox"/> Pull/Steering Bar <b>*cannot be used with rear steering</b>

<input type="checkbox"/> <b>AM-12Small HAND-FOOT TRYKE</b> (with 3" crank arms & 12" wheels) <b>50-HFC-0110</b> Arm length 12-14", leg length 16-20" <b>Now comes standard with rear steering!</b>
<input type="checkbox"/> Option #1 – Standard Seating System: Bench Seat & 1600 Simple Seat Back
<input type="checkbox"/> Option #2 – Upgraded Seating System: Snappy Seat System Snappy Seat Accessories: <input type="checkbox"/> Notched Seat Bottom Cushion* <input type="checkbox"/> Notched Wide Bottom Cushion* <input type="checkbox"/> XL Seat Back <input type="checkbox"/> Laterals <input type="checkbox"/> H-Harness: <input type="checkbox"/> 8" <input type="checkbox"/> 11.5"
<input type="checkbox"/> Option #3 – Upgraded Seating System: Blue Bucket Seat
AM-12S Accessories: <input type="checkbox"/> Pull/Steering Bar <input type="checkbox"/> Toe Pedal Leveler Pulley <input type="checkbox"/> Exercise pedals <input type="checkbox"/> Push Bar <input type="checkbox"/> Heavy-Duty Push Bar <b>*cannot be used with rear steering</b>

<input type="checkbox"/> <b>AM-12 HAND-FOOT TRYKE</b> (with 4" crank arms & 12" wheels) <b>50-HFC-0210</b> Arm length 14-23", leg length 21-29" <b>Now comes standard with rear steering!</b>
<input type="checkbox"/> Option #1 – Standard Seating System: Saddle Seat/1600 Simple Seat Back Seat Bottom Alternates: <input type="checkbox"/> Bench Seat <input type="checkbox"/> Pommel Saddle Seat <input type="checkbox"/> Medium (11.5x10") <input type="checkbox"/> Large(14x13") Trunk Support Upgrades: <input type="checkbox"/> ProSeries Seatback ↴ <input type="checkbox"/> ProSeries Head Rest <input type="checkbox"/> ProSeries Full Padded Back (includes push grip, backpad & 2 laterals) <input type="checkbox"/> ProSeries Lumbar Pad <input type="checkbox"/> ProSeries Back Spacer Kit <input type="checkbox"/> H-Harness: <input type="checkbox"/> 8" <input type="checkbox"/> 11.5" <input type="checkbox"/> AM Full Padded Back
<input type="checkbox"/> Option #2 – Upgraded Seating System: Gray Bucket Seat Bucket Seat Accessories: <input type="checkbox"/> H-Harness: <input type="checkbox"/> 8" <input type="checkbox"/> 11.5"
<input type="checkbox"/> Option #3 – Upgraded Seating System: Snappy Seat System Snappy Seat Accessories: <input type="checkbox"/> Notched Seat Bottom Cushion <input type="checkbox"/> Notched Wide Bottom Cushion <input type="checkbox"/> XL Seat Back <input type="checkbox"/> Laterals <input type="checkbox"/> H-Harness: <input type="checkbox"/> 8" <input type="checkbox"/> 11.5"
AM-12 Accessories: <input type="checkbox"/> Pull/Steering Bar <input type="checkbox"/> Toe Pedal Leveler Pulley <input type="checkbox"/> Exercise pedals <input type="checkbox"/> Knee Separator: <input type="checkbox"/> 2" <input type="checkbox"/> 5"

Recipient's Name:

**AM-16 HAND-FOOT TRYKE** (with 16" wheels) **50-HFC-0411** Arm length 18-27", leg length 24-36"

**Now comes standard with rear steering!**

Option #1 – Standard Seating System: Saddle Seat/Gray Plastic Back

Seat Bottom Upgrades:  Bench Seat  Padded Tractor Seat (used with ProSeries Seatback only)

Pommel Saddle Seat  Medium (11.5x10")  Large(14x13")

Trunk Support Upgrades:  ProSeries Seatback ↴  ProSeries Head Rest  ProSeries Full Padded Back  
(includes push grip, backpad & 2 laterals)

ProSeries Lumbar Pad  ProSeries Back Spacer Kit

H-Harness:  8"  11.5"  AM Full Padded Back

Option #2 – Standard Seating System: Saddle Seat/1600 Simple Seat Back

Seat Bottom Upgrades:  Bench Seat  Padded Tractor Seat (used with ProSeries Seatback only)

Pommel Saddle Seat  Medium (11.5x10")  Large(14x13")

Trunk Support Upgrades:  ProSeries Seatback ↴  ProSeries Head Rest  ProSeries Full Padded Back  
(includes push grip, backpad & 2 laterals)

ProSeries Lumbar Pad  ProSeries Back Spacer Kit

H-Harness:  8"  11.5"  AM Full Padded Back

Option #3 – Upgraded Seating System: Gray Bucket Seat

Bucket Seat Accessories:  H-Harness:  8"  11.5"

AM-16 Accessories:  Pull/Steering Bar  Toe Pedal Leveler Pulley  Exercise pedals  XL Exercise pedals  
 Knee Separator:  2"  5"  Knee Separator 9" Extender Tube

**New! AM-20 HAND-FOOT TRYKE** (Recumbent Single Speed) – **50-HFC-0610** Arm length 18-30", leg length 32-46"

AM-20 Accessories:  Exercise Pedals  XL Exercise Pedals  Expanding Pedals

Wheelchair Seat Bar Ends  Swing Away Arms

Therapist Assembly Notes or Comments (if any):

Recipient's Name: \_\_\_\_\_

## HAND TRYKES

<input type="checkbox"/> <b>New! AM-10 HAND TRYKE</b> – (4" Crank Arms & Permanent Foot Platform) <b>50-HC-0105</b> Arm Length 12-16"
<input type="checkbox"/> Option #1 – Standard Seating System: Blue Bucket Seat
<input type="checkbox"/> Option #2 – Upgraded Seating System: Snappy Seat System Snappy Seat Accessories: <input type="checkbox"/> Notched Seat Bottom Cushion <input type="checkbox"/> Notched Wide Bottom Cushion <input type="checkbox"/> XL Seat Back <input type="checkbox"/> Laterals <input type="checkbox"/> H-Harness: <input type="checkbox"/> 8" <input type="checkbox"/> 11.5"
AM-10 Accessories: <input type="checkbox"/> Push Bar <input type="checkbox"/> Heavy-Duty Push Bar <input type="checkbox"/> Pull/Steering Bar <input type="checkbox"/> Rear Steering Kit

<input type="checkbox"/> <b>AM-12Small HAND TRYKE</b> – (12" Front Wheel & Permanent Foot Platform) <b>50-HC-0210S</b> Arm Length 12-16" <b>Now comes standard with rear steering!</b>
<input type="checkbox"/> Option #1 – Standard Seating System: Bench Seat & 1600 Simple Seat Back Trunk Support Upgrades: <input type="checkbox"/> AM Full Padded Back <input type="checkbox"/> H-Harness: <input type="checkbox"/> 8" <input type="checkbox"/> 11.5" <input type="checkbox"/> Back Spacer Kit
<input type="checkbox"/> Option #2 – Upgraded Seating System: Snappy Seat System Snappy Seat Accessories: <input type="checkbox"/> Notched Seat Bottom Cushion* <input type="checkbox"/> Notched Wide Bottom Cushion* <input type="checkbox"/> XL Seat Back <input type="checkbox"/> Laterals <input type="checkbox"/> H-Harness: <input type="checkbox"/> 8" <input type="checkbox"/> 11.5"
<input type="checkbox"/> Option #3 – Upgraded Seating System: Blue Bucket Seat
AM-12Small Accessories: <input type="checkbox"/> Pull/Steering Bar <input type="checkbox"/> Push Bar <input type="checkbox"/> Heavy-Duty Push Bar <input type="checkbox"/> Knee Separator: <input type="checkbox"/> 2" <input type="checkbox"/> 5" <b>*cannot be used with rear steering</b>

<input type="checkbox"/> <b>AM-12 HAND TRYKE</b> – (12" Front Wheel & Permanent Foot Platform) <b>50-HC-0210</b> Arm Length 14-23" <b>Now comes standard with rear steering!</b>
<input type="checkbox"/> Option #1 – Standard Seating System: Saddle Seat/1600 Simple Seat Back Seat Bottom Alternates: <input type="checkbox"/> Bench Seat <input type="checkbox"/> Pommel Saddle Seat <input type="checkbox"/> Medium (11.5x10") <input type="checkbox"/> Large(14x13") Trunk Support Upgrades: <input type="checkbox"/> ProSeries Seatback ↴ <input type="checkbox"/> ProSeries Head Rest <input type="checkbox"/> ProSeries Full Padded Back (includes push grip, backpad & 2 laterals) <input type="checkbox"/> ProSeries Lumbar Pad <input type="checkbox"/> ProSeries Back Spacer Kit <input type="checkbox"/> H-Harness: <input type="checkbox"/> 8" <input type="checkbox"/> 11.5" <input type="checkbox"/> AM Full Padded Back
<input type="checkbox"/> Option #2 – Upgraded Seating System: Gray Bucket Seat Bucket Seat Accessories: <input type="checkbox"/> H-Harness: <input type="checkbox"/> 8" <input type="checkbox"/> 11.5"
<input type="checkbox"/> Option #3 – Upgraded Seating System: Snappy Seat System Snappy Seat Accessories: <input type="checkbox"/> Notched Seat Bottom Cushion* <input type="checkbox"/> Notched Wide Bottom Cushion* <input type="checkbox"/> XL Seat Back <input type="checkbox"/> Laterals <input type="checkbox"/> H-Harness: <input type="checkbox"/> 8" <input type="checkbox"/> 11.5"
AM-12 Accessories: <input type="checkbox"/> Pull/Steering Bar <input type="checkbox"/> Knee Separator: <input type="checkbox"/> 2" <input type="checkbox"/> 5" <b>*cannot be used with rear steering</b>

Recipient's Name:

**AM-16 HAND TRYKE** (16" wheels) **50-HFC-0411** Arm length 18-27"

**Now comes standard with rear steering!**

Option #1 – Standard Seating System: Saddle Seat/Gray Plastic Back

- Seat Bottom Upgrades:  Bench Seat  Padded Tractor Seat (used with ProSeries Seatback only)  
 Pommel Saddle Seat  Medium (11.5x10")  Large(14x13")
- Trunk Support Upgrades:  ProSeries Seatback ↴  ProSeries Head Rest  ProSeries Full Padded Back  
(includes push grip, backpad & 2 laterals)  
 ProSeries Lumbar Pad  ProSeries Back Spacer Kit  
 H-Harness:  8"  11.5"  AM Full Padded Back

Option #2 – Standard Seating System: Saddle Seat/1600 Simple Seat Back

- Seat Bottom Upgrades:  Bench Seat  Padded Tractor Seat (used with ProSeries Seatback only)  
 Pommel Saddle Seat  Medium (11.5x10")  Large(14x13")
- Trunk Support Upgrades:  ProSeries Seatback ↴  ProSeries Head Rest  ProSeries Full Padded Back  
(includes push grip, backpad & 2 laterals)  
 ProSeries Lumbar Pad  ProSeries Back Spacer Kit  
 H-Harness:  8"  11.5"  AM Full Padded Back

Option #3 – Upgraded Seating System: Gray Bucket Seat

Bucket Seat Accessories:  H-Harness:  8"  11.5"

AM-16 Accessories:  Pull/Steering Bar  
 Knee Separator:  2"  5"  Knee Separator 9" Extender Tube

**COMMUNITY CRUISER HAND TRYKE** – (24" Front Wheel & 3-Speed Drive Train) **50-HC-1424** Arm length 22-26"

CC Accessories:  Wheelchair Seat Bar Ends  Swing Away Arms  Remote location break and shifting kit

**HP-1000 RECREATION HAND TRYKE** – (7-Speed Drive Train) **50-HC-1000** Arm length 22-26"


HP-1000 Accessories:  Wheelchair Seat Bar Ends  Swing Away Arms

Therapist Assembly Notes or Comments (if any):

# Amtryke Sizing Chart

TRYKE TYPE (How will the tryke be propelled?)	RIDER LEG LENGTH (Inches from center of hip to bottom of shoe.)	RIDER ARM LENGTH (Inches from middle of shoulder to center of digit crease.)	MODEL	RIDER WEIGHT (Pounds)	RIDER MAX HEIGHT (Inches)	TRYKE WEIGHT (Pounds)	WHEEL SIZE (Inches)	TTRYKE HEIGHT (Inches)	TRYKE LEGNTH (Inches)	TRYKE WIDTH (Inches)
<b>Hand &amp; Foot</b>	14-21	12-16	AM-10	175	40	45	11	24	38	21
	18-22	12-16	AM-12S	150	40	45	12	27	38	24
	21-29	14-23	AM-12	150	47	45	12	36	60	32
	24-36	18-27	AM-16	175	66	55	16	36	68	33
	32-46	18-30	AM-20	250	74	80	20	45	72	32
<b>Foot</b>	14-21	7-20	1410	175	40	45	11	24	38	21
	17-22	9-13	Snappy	150	42	45	12	30	38	24
	21-26	14-22	1412	125	42	72	12	13	43	27
	24-32	14-22	1416/1616	175	60	74	16	49	58	30
	27-35	18-28	1420/1620	250	68	74	20	50	64	30
	28-43	20-32	1420XL	275	76	89	20	43	72	29
	28-36	22-30	2722	275	72	76	24	42	60	30
	30-41	20-28	JT-2000/ JT2300USS	250	74	80	20	48	72	32
<b>Hand</b>	22-41	22-26	1424	250	72	74	20	45	75	32
	N/A	22-26	HP-1000	250	74	80	Front 16 Rear 20	45	72	32
<b>All trykes in the Hand &amp; Foot section can be converted to Hand trykes – except the AM-20.</b>										

- A Center of Shoulder
- B Center of Elbow
- C Center of Digit Crease
- D Center of Hip
- E Center of Knee
- F Bottom of Foot



RIDER'S MEASUREMENTS	
<b>Arm Measurements (inches) Total Length</b>	
Left	A to B + B to C = _____
Right	A to B + B to C = _____
Trunk	A to D = _____
<b>Leg Measurements (inches) Total Length</b>	
Left	D to E + E to F = _____
Right	D to E + E to F = _____

HELMET SIZING	
<b>Youth Sizes</b>	<b>Head Circumference Inches</b>
XXS	18.5 to 19.5
XS	20 to 21
S/M	21.5 to 22.5
L/XL	22.75 to 24.5
<b>Adult Sizes</b>	<b>Head Circumference Inches</b>
S/M	22 to 23 5/8
L/XL	23 5/8 to 25 3/4